



PRODUCT DATA SHEET

POLYVIEW[®] IHC reagent (rabbit-HRP)

ADI-950-115

Highest sensitivity, lowest background nanopolymer detection reagent for use with HIGHDEF[®] chromogens.

Product Number/Sizes

ADI-950-115-0100	100 tests
ADI-950-115-0025	25 tests

- Biotin-free nanopolymer detection circumvents endogenous biotin background
- High intensity color development delivers sharp, crisp staining
- Customize with HIGHDEF[®] chromogens of your choice (sold separately)

The POLYVIEW[®] IHC reagent (rabbit-HRP) is a non-biotin one-step detection reagent suitable for demonstrating antigens in formalin-fixed paraffin-embedded tissues and frozen sections. The detection reagent may also be used with blood smears, cytosmears, and cell preparations.

Product Specifications

APPLICATIONS:	IHC
KIT/SET CONTAINS:	Anti-rabbit HRP detection reagent
SHIPPING:	Blue Ice Not Frozen
LONG TERM STORAGE:	+4°C
TECHNICAL INFO/PRODUCT NOTES:	

The POLYVIEW[®] IHC reagent (rabbit-HRP) has been developed by directly labeling anti-rabbit immunoglobulins with enzymes using a proprietary tandem hyperlabelling technology. This ensures consistent and reproducible immunodetection of rabbit antibodies against nuclear, cytoplasmic and membrane antigens in different types of tissues. The single step methodology enables faster staining procedures than traditional two-step methods using biotin and avidin/streptavidin conjugates, with significantly lower background. The POLYVIEW[®] IHC reagent (rabbit-HRP) is suitable for use with rabbit antibodies, both monoclonal and polyclonal. The reagents can be used for manual staining or with automated staining instruments and are well suited for multiplex immunohistochemical staining assays.

Application Note:

Compatibility of POLYVIEW[®] IHC Detection Reagents and HIGHDEF[®] IHC Chromogens with Automated IHC Slide Stainers

全国代理




www.boppard.cn
info@boppard.cn

北京 Tel: 010 85804838
上海 Tel: 021 62884751
广州 Tel: 020 87326381
香港 Tel: 852 27999019


For Research Use Only, Not for Human

GLOBAL HEADQUARTERS

Enzo Life Sciences, Inc.
10 Executive Blvd
Farmingdale, NY 11735
USA
T 1-800-942-0430
T 1-631-694-7070
F 1-610-941-9252
E info-usa@enzolifesciences.com
www.enzolifesciences.com



Expression of Ki67 in mouse intestine. Immunohistochemistry staining of Ki67 (A, B, C) in formalin-fixed paraffin-embedded mouse intestine using POLYVIEW[®] IHC reagent (rabbit-HRP) (Prod. No. ADI-950-115) and HIGHDEF[®] IHC chromogen substrate (DAB, stable) (Prod. No. ADI-950-212) from Enzo Life Sciences (A), Dako (B), and Vector Labs (C).



Expression of p53 in mouse intestine. Immunohistochemistry staining of p53 (A, B, C) in formalin-fixed paraffin-embedded mouse intestine using POLYVIEW[®] IHC reagent (rabbit-HRP) (Prod. No. ADI-950-115) and HIGHDEF[®] IHC chromogen substrate (DAB, stable) (Prod. No. ADI-950-212) from Enzo Life Sciences (A), Dako (B), and Vector Labs (C).

Revised 31-Mar-16

全国代理



宝柏 · 中国

www.boppard.cn

info@boppard.cn

北京 Tel: 010 85804838

上海 Tel: 021 62884751

广州 Tel: 020 87326381

香港 Tel: 852 27999019

For Research Use Only, Not for Human

GLOBAL HEADQUARTERS

Enzo Life Sciences, Inc.
10 Executive Blvd
Farmingdale, NY 11735
USA

T 1-800-942-0430

T 1-631-694-7070

F 1-610-941-9252

E info-usa@enzolifesciences.com

www.enzolifesciences.com