<For Research Use Only>

Wakofor Genetic Research
PCR Purification Kit Wako**[Introduction]**

PCR Purification Kit *Wako* is used for removal of primers, primer dimers, dNTPs, enzymes such as DNA polymerase and restriction enzyme from reaction buffer of various molecular biology experiments.

[Features]

- 1) Purification of the DNA fragment in about 10 minutes
- 2) Purification of 80 bp–30 kbp DNA
- 3) Removal of primers, primer dimers, dNTPs, and enzymes

[Kit contents] (30 purifications)

	Contents
(1) Binding Buffer	5 mL × 1 vial
(2) Washing Buffer ^{*1}	12 mL × 1 vial
(3) Elution Buffer	2 mL × 1 tube
(4) Spin Column	30 columns × 1 pack
(5) 2 mL Collection Tube	30 tubes × 1 pack
(6) 1.5 mL Collection Tube ^{*2}	30 tubes × 1 pack

※ 1. Add 18 mL of 99.5% ethanol into the vial of Washing Buffer before use of this kit.

※ 2. Thirty tubes of 1.5 mL Collection Tube are included in this kit for the DNA elution step of [Procedure] (1) Steps 8. ~ 9. In the case of [Procedure] (2), you have to prepare 30 new tubes for use of Step 1.

[Additional materials required]

Equipments :

- 1) Microcentrifugation tube
- 2) PCR centrifugation tube
- 3) Micropipette
- 4) Pipette tips
- 5) UV Gel illuminator

Reagents :

- 1) 99.5% Ethanol

[Precautions]

1. Equipments

PCR tube, microcentrifugation tube, and pipette tips should be sterilized by autoclaving. There is no problem with using commercially available DNase and RNase free products. We recommend that you use gloves and a mask to avoid RNase contamination.

2. Reagents

Additional reagents excluding the kit are sterilized by autoclaving or filtration.

※ 1. Add 18 mL of 99.5% ethanol into the vial of Washing Buffer before use of this kit.

[Procedures](1) DNA purification from PCR solution ($\leq 50 \mu\text{L}$)

1. Add 130 μL of Binding Buffer into the PCR microcentrifugation tube including DNA solution after PCR and mix by pipetting gently.
2. Set Spin Column in a 2 mL Collection Tube.
3. Add the total solution of step 1 into the Spin Column.
4. Incubate for 1 minute at room temperature and centrifuge at 10,000 rpm for 30 seconds at room temperature.
5. Add 700 μL of Washing Buffer^{*1} into the Spin Column and centrifuge at 10,000 rpm for 30 seconds at room temperature.
6. Remove Spin Column^{*3} from the 2 mL Collection Tube and dispose the solution in the 2 mL Collection Tube.
7. Reset the Spin Column in the 2 mL Collection Tube and centrifuge at 12,000 ~ 14,000 rpm for 3 minutes at room temperature.^{*4}
8. Transfer the Spin Column into a 1.5 mL Collection Tube and add 20 μL of Elution Buffer.^{*5}
9. Incubate for 3 minutes^{*5} at room temperature and centrifuge at 10,000 rpm for 1 minute at room temperature.
10. Use the purified DNA for the following applications.
nuclease digestion, southern blot analysis, labeling, cloning, sequencing, *in vitro* transcription, etc.

※ 3. Keep the Spin Column in this step.

※ 4. This procedure is to remove the remaining alcohol solution from the Spin Column.

※ 5. In the case of low DNA recovery, we recommend increasing Elution Buffer to 40 ~ 80 μL , or to extending incubation time from 3 minutes to 10 minutes.


(2) DNA purification from PCR solution (50 ~ 200 μL)

1. Transfer DNA solution after PCR into a new 1.5 mL microcentrifugation tube.
2. Add 300 μL of Binding Buffer and mix by pipetting gently.
3. Set the Spin Column in a 2 mL Collection Tube.
4. Add the total solution of step 1 into the Spin Column.
5. Incubate for 1 minute at room temperature and centrifuge at 10,000 rpm for 30 seconds at room temperature.
6. Add 700 μL of Washing Buffer^{*1} and centrifuge at 10,000 rpm for 30 seconds at room temperature.
7. Remove the Spin Column^{*6} from the 2 mL Collection Tube and dispose the solution in the 2 mL Collection Tube.
8. Reset the Spin Column in the 2 mL Collection Tube and centrifuge at 12,000 ~ 14,000 rpm for 3 minutes at room temperature.^{*7}
9. Transfer the Spin Column into a 1.5 mL Collection Tube and add 30 μL of Elution Buffer into the Spin Column.^{*8}
10. Incubate for 3 minutes^{*8} at room temperature and centrifuge at 10,000 rpm for 1 minute at room temperature.
11. Use the purified DNA for following applications.
nuclease digestion, southern blot analysis, labeling, cloning, sequencing, *in vitro* transcription, etc.

※ 6. Please keep the Spin Column in this step.

※ 7. This procedure is to remove the remaining alcohol solution from the Spin Column.

※ 8. In the case of low DNA recovery, we recommend increasing the Elution Buffer to 60 ~ 100 μL , or extending the incubation time from 3 minutes to 10 minutes.

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