<For Research Use Only>

Code No. 298-68001 (10 reactions)

Wako

for genetic research
Target mRNA Cloning Kit Wako


[Introduction]

Target mRNA Cloning Kit *Wako* is used for cDNA amplification from a small amount of mRNA. This product is used for the cloning of target mRNA bound with RISC (RNA-induced silencing complex) using an Ago IP RNA fraction.

[Features]

- 1) Independent cDNA amplification against the length of mRNA
- 2) High efficiency of cDNA amplification
- 3) Investigation of microRNA targeted mRNA

[Outline of kit procedures]

[Kit contents] (10 Reactions)

Target mRNA Cloning Kit *Wako* includes a total of 14 components

	Content
(1) dT(20)-RT Primer (20 pmol/ μ L)	10 μ L \times 1 tube
5'-GACCATATGACGAGATCCGAGCTTCTTTTTTTTTTTTTTTTTTTT-3' (45mer)	
(2) 10 \times RT Buffer	20 μ L \times 1 tube
(3) dNTP Mixture Solution (2.5 mmol/L each)	20 μ L \times 1 tube
(4) RNase Inhibitor (20 U/ μ L)	10 μ L \times 1 tube
(5) Reverse Transcriptase (200 U/ μ L)	10 μ L \times 1 tube
(6) 0.5 mol/L EDTA Solution	20 μ L \times 1 tube
(7) 1 mol/L Tris-HCl Solution	200 μ L \times 1 tube
(8) Ethachinmate	30 μ L \times 1 tube
(9) 10 mol/L Ammonium Acetate Solution	480 μ L \times 1 tube
(10) 2nd Strand Synthesis Primer (20 pmol/ μ L)	10 μ L \times 1 tube
5'-CATGGTATCGACGAGACTGAGGCTGNNNNNNNNB-3' (35mer)	
(11) 5' PCR Primer (20 pmol/ μ L)	10 μ L \times 1 tube
5'-CATGGTATCGACGAGACTGAG-3' (21mer)	
(12) 3' PCR Primer (20 pmol/ μ L)	10 μ L \times 1 tube
5'-GACCATATGACGAGATCCGAG-3' (21mer)	
(13) 5' Colony PCR Primer (5 pmol/ μ L)	960 μ L \times 1 tube
5'-GTATCGACGAGACTGAGGCTG-3' (21mer)	
(14) 3' Colony PCR Primer (5 pmol/ μ L)	960 μ L \times 1 tube
5'-CGAGATCCGAGCTTCTTTTT-3' (21mer)	

[Additional materials required]

Reagents :

[1] Thermostable DNA polymerase for PCR

We recommend HotGoldstar™ DNA polymerase (Code No. 314-80251, 310-80253) for suppression of non-specific amplification. HotGoldstar™ DNA polymerase is used for hot start PCR and activated by 10 minutes incubation at 95°C before PCR cycles. DNA fragments, which were amplified by HotGoldstar™ DNA polymerase, are 3' protruding end of adenine base. These fragments therefore are able to directly insert into T-vectors.

[2] PCR Purification Kit *Wako* (Wako Cat. # 298-67901)

- 3) Sterile distilled water
- 4) 70%, 99.5% Ethanol
- 5) 0.2 mol/L Sodium hydroxide solution
- 6) Cloning vector
- 7) *E. coli* competent cells
- 8) Ligation reagents
- 9) PBS (-)
- 10) Phenol : Chloroform : Isoamyl alcohol (25:24:1)
- 11) Chloroform
- 12) Trypsin-EDTA Solution (for adherent cells)

Equipment :

- 1) 1.5 mL Micro centrifugation tube
- 2) Micro pipette
- 3) Pipette tip
- 4) Tubes rotator
- 5) High speed microcentrifuge (Max > 14,500 rpm, 20,000 \times g)
- 6) Vortex mixer

[Important information before experiment]

1. Equipment
PCR tube, microcentrifugation tube, pipette tip must be sterilized by autoclave. There is no problem with using commercially available RNase free products. We recommend that you must use gloves and a mask to avoid contamination of RNase.
2. Reagents
Additional reagents excluding kit must be sterilized by autoclave or filtration.

[Procedure]

«Preparation of Ago2 IP RNA»

The protocol of this step is same with following products until the step of 4.-b) of this instruction manual.

- microRNA Isolation Kit, Human Ago2 (Wako cat.# 292-66701)
- microRNA Isolation Kit, Mouse Ago2 (Wako cat.# 292-67301)

This kit is covering the following step of «Reverse Transcription».

1. Preparation of cell pellet

For adherent cells

- 1) Culture the cells ($5 \times 10^6 \sim 1 \times 10^7$ cells).^{※1}
- 2) Remove the medium from culture dishes.
- 3) Wash the cells twice with PBS (-).
- 4) Collect the cells using a scraper or Trypsin-EDTA solution from dishes. In the case of using Trypsin-EDTA solution, add medium including serum after cell collection.
- 5) Centrifuge cell suspension at 200 \times g for 5 minutes at 4°C.
- 6) Remove the supernatant.
- 7) Suspend cell pellet in 1 mL of PBS (-).
- 8) Transfer all of suspension into a new 1.5 mL micro centrifugation tube.
- 9) Repeat steps 5) to 7).
- 10) Repeat steps 5) to 6).

For suspension cells

- 1) Culture the cells ($5 \times 10^6 \sim 1 \times 10^7$ cells).^{※1}
- 2) Centrifuge cell suspension at $200 \times g$ for 5 minutes at 4°C.
- 3) Remove the supernatant.
- 4) Suspend cell pellet in 1 mL of PBS (-).
- 5) Transfer all of suspension into a new 1.5 mL micro centrifugation tube.
- 6) Repeat steps 2) to 4).
- 7) Repeat steps 2) to 3).

※1. Maximum cell numbers are 1×10^7 . Non-specific immunoprecipitation and the decrease of microRNA purification efficiency will be caused by the use of more than $1 \times 10^7 \times$ cells.

2. Prewashing of Anti Human/Mouse Ago2 Antibody Beads

- 1) Transfer 50 μ L of Anti Human/Mouse Ago2 Antibody Beads Solution into a new 1.5 mL micro centrifugation tube.
- 2) Centrifuge beads suspension at $3,000 \times g$ for 30 seconds at 4°C.
- 3) Remove the supernatant.
- 4) Suspend the beads pellet by using 1 mL of Cell Lysis Solution.
- 5) Repeat steps 2) and 3). (Prewashed Beads)

3. Preparation of Cell Lysate

- 1) Add 1 mL of Cell Lysis Solution into a tube of cell pellet and suspend cells by pipetting.
- 2) Incubate for 10 minutes on ice.
- 3) Centrifuge the cell lysate at $20,000 \times g$ for 20 minutes at 4°C.
- 4) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. (Cell Lysate)

4. Purification of microRNA fraction

a) Antigen-Antibody Reaction (Immunoprecipitation)

- 1) Add Cell Lysate into a tube of Prewashed Beads.
- 2) Mix by vortex mixer.
- 3) Mix by rotator for minimum 2 hours at 4°C.
- 4) Centrifuge at $3,000 \times g$ for 30 seconds at 4°C.
- 5) Remove the supernatant. (Reacted Beads)

b) Washing after Antigen-Antibody Reaction and elution of Ago2 protein

- 1) Add 1 mL of Cell Lysis Solution into a tube of Reacted Beads.
- 2) Mix by vortex mixer.
- 3) Centrifuge at $3,000 \times g$ for 30 seconds at 4°C.
- 4) Remove the supernatant.
- 5) Repeat steps 1) to 4) at twice. (Washed Beads)
- 6) Add 50 μ L of Elution Solution^{*1} into a tube of Washed Beads and mix by vortex mixer.
- 7) Centrifuge at $3,000 \times g$ for 30 seconds at 4°C.
- 8) Transfer the supernatant^{*2} into a new 1.5 mL micro centrifugation tube. (Protein Solution)

c) Purification of microRNA

- 1) Add 1 mL of ISOGEN^{*3} reagent into a tube of Protein Solution and mix by vortex mixer.
- 2) Add 300 μ L of chloroform and mix by vortex mixer for 1 minute. Incubate for 30 minutes on ice.
- 3) Centrifuge at $18,800 \times g$ for 10 minutes at 4°C.
- 4) Transfer the supernatant^{*4} into a new 2.0 mL micro centrifugation tube. Add 1 μ L of Ethachinmate and equal volume of Isopropanol. Mix by vortex mixer.
- 5) Centrifuge at $18,800 \times g$ for 10 minutes at 4°C.
- 6) Remove the supernatant.
- 7) Add 500 μ L of 70% Ethanol and mix by vortex mixer.
- 8) Centrifuge at $18,800 \times g$ for 10 minutes at 4°C.
- 9) Repeat steps 6) to 8).

10) Remove the supernatant.

- 11) Open the cap of the 2.0 mL micro centrifugation tube and dry the pellet for 20 minutes at room temperature.
- 12) Dissolve the pellet in 13 μ L of sterile distilled water.^{*5} (Ago2 IP RNA Solution)

- *1. Components of the Elution Solution will be precipitated during storage at 4°C. In this case, Elution Solution should be warmed to room temperature and mixed by a vortex mixer. We recommend confirmation of whether precipitation is completely dissolved or not.
- *2. To avoid the contamination of non-specific nucleic acids, carefully pipette the supernatant without absorption of beads pellet.
- *3. This reagent is available from NIPPON GENE CO., LTD.. (Cat No. 311-02501) This reagent mainly is included guanidine thiocyanate and phenol.
- *4. Carefully pipette only the supernatant.
- *5. Store the sample solution at -80°C.

◀Reverse Transcription▶

The kit components of this product are used for the following protocol. We recommend the use of PCR Purification Kit *Wako* (Wako cat.# 298-67901), which is sold separately, for the steps of <Reverse Transcription>, <2nd strand cDNA synthesis> and <PCR Amplification>.

- 1) Transfer 13 μ L of Ago2 IP RNA Solution into a PCR tube and add 1 μ L of dT(20)-RT Primer.
- 2) Incubate for 3 minutes at 70°C and for 2 minutes on ice.
- 3) Mix the following reagents in a PCR tube on ice.

<u>Ago2 IP RNA Solution</u> + dT(20)-RT Primer	14 μ L
10 \times RT Buffer	2 μ L
dNTP Mixture Solution	2 μ L
RNase Inhibitor	1 μ L
Reverse Transcriptase	1 μ L
	20 μ L

- 4) Incubate for 10 minutes at 42°C. (Reverse transcription)
- 5) Add 2 μ L of 0.5 mol/L EDTA and mix gently on ice. (Stop reaction at RT)
- 6) Add 8 μ L of 0.2 mol/L NaOH (50 μ L of 2 mol/L NaOH + 450 μ L of sterile distilled water) (Not included in this kit).
- 7) Mix and incubate for 30 minutes at 65°C. (Alkaline hydrolysis of template mRNA)
- 8) Incubate for 1 minute on ice.
- 9) Add 20 μ L of 1 mol/L Tris-HCl Solution and mix.
- 10) Use the PCR Purification Kit *Wako* for purification of nucleic acid by the following procedure.

◀Procedure of the PCR Purification Kit *Wako* (Wako cat.# 298-67901)▶

The kit components of PCR Purification Kit *Wako* : 「※」

- 1) Add 130 μ L of Binding Buffer[※] and mix.
- 2) Set Spin Column[※] in a 2 mL Collection Tube[※].
- 3) Transfer the total solution of step 1) into the Spin Column[※].
- 4) Incubate for 1 minute at room temperature.
- 5) Centrifuge at $9,600 \times g$ for 1 minute at room temperature.
- 6) Remove the Spin Column[※] from the 2 mL Collection Tube and dispose of the solution in the 2 mL Collection Tube.
- 7) Add 700 μ L of Washing Buffer[※] into the Spin Column[※].

- 8) Repeat steps 5) to 6).
- 9) Reset the Spin Column[®] in the 2 mL Collection Tube[®] and centrifuge at 18,800 × g for 3 minutes at room temperature.
- 10) Transfer the Spin Column[®] into a 1.5 mL Collection Tube[®].
- 11) Add 50 μL of Elution Buffer[®] into the Spin Column[®].
- 12) Incubate for 3 minutes at room temperature.
- 13) Centrifuge at 9,600 × g for 1 minute at room temperature.
- 14) Add 1 μL of Ethachinmate and mix.
- 15) Add 12 μL of 10 mol/L Ammonium Acetate Solution and mix.
- 16) Add 125 μL of 99.5% Ethanol and mix.
- 17) Centrifuge at 18,800 × g for 10 minutes at 4 °C.
- 18) Remove the supernatant.
- 19) Add 500 μL of 70% Ethanol and mix by vortex mixer.
- 20) Centrifuge at 18,800 × g for 10 minutes at 4 °C.
- 21) Repeat steps 18) to 20).
- 22) Remove the supernatant.
- 23) Dry the pelle.
- 24) Dissolve the pellet in 35.5 μL of sterile distilled water. (1st Strand cDNA Solution)

«2nd strand cDNA synthesis»

We recommend the use of HOTGoldstar[™] DNA Polymerase (EUROGENTEC), which is sold separately.

- 1) Transfer 35.5 μL of 1st Strand cDNA Solution into a PCR tube and mix the following reagents in the PCR tube on ice.

1st Strand cDNA Solution	35.5 μL
2nd Strand Synthesis Primer	1 μL
10 × PCR Reaction Buffer (EUROGENTEC)	5 μL
2.5 mmol/L dNTP Mixture (EUROGENTEC)	4 μL
25 mmol/L MgCl ₂ (EUROGENTEC)	4 μL
HOTGoldstar [™] DNA Polymerase (EUROGENTEC)	0.5 μL
Total	50 μL

- 2) Incubate at following condition. (2nd strand cDNA synthesis)

95 °C	10 minutes
30 °C	1 minute
40 °C	30 seconds
50 °C	30 seconds
60 °C	30 seconds
72 °C	2 minutes
95 °C	1 minute
氷上	

- 3) Use the PCR Purification Kit *Wako* (Protocol steps 1)–22)) for purification of nucleic acid.
- 4) Dry the pellet.
- 5) Dissolve the pellet in 34.5 μL of sterile distilled water. (2nd Strand cDNA Solution)

«PCR Amplification»

- 1) Transfer 34.5 μL of 2nd Strand cDNA Solution into a PCR tube and mix the following reagents in the PCR tube on ice.

<u>2nd Strand cDNA Solution</u>	34.5 μL
5' PCR Primer	1 μL
3' PCR Primer	1 μL
10 × PCR Reaction Buffer (EUROGENTEC)	5 μL
2.5 mmol/L dNTP Mixture (EUROGENTEC)	4 μL
25 mmol/L MgCl ₂ (EUROGENTEC)	4 μL
HOTGoldstar [™] DNA Polymerase (EUROGENTEC)	0.5 μL
Total	50 μL

- 2) Incubate at the following conditions. (PCR amplification)

95 °C	10 minutes	
↓		
95 °C	20 seconds	} 24–26 cycles
60 °C	10 seconds	
72 °C	30 seconds	
↓		
72 °C	2 minutes	

- 3) Incubate for 1 minute on ice.
- 4) Use the PCR Purification Kit *Wako* (Protocol steps 1)–22)) for purification of nucleic acid.
- 5) Dry the pellet.
- 6) Dissolve the pellet in 4 μL of sterile distilled water. (cDNA Solution)

«Transformation of *E. coli*»

- 1) Ligate cDNA Solution into T-vector^{*5}.
- 2) Transform the constructed plasmids into *E. coli* competent cells. Select transformants by using LB agar medium including antibiotics.

*5. We recommend checking the DNA length and concentration of cDNA Solution before TA cloning by using Bioanalyzer.

«Insert check by colony PCR»

- 1) Mix the following reagents in a PCR tube on ice.

5' Colony PCR Primer	1 μL
3' Colony PCR Primer	1 μL
10 × PCR Reaction Buffer (EUROGENTEC)	1 μL
2.5 mmol/L dNTP Mixture (EUROGENTEC)	1 μL
25 mmol/L MgCl ₂ (EUROGENTEC)	0.8 μL
Sterile distilled water	5.1 μL
HOTGoldstar [™] DNA Polymerase (EUROGENTEC)	0.1 μL
Total	10 μL

- 2) Suspend the colony in this reaction solution, and incubate at the following conditions.

95 °C	10 minutes	
↓		
95 °C	20 seconds	} 30 cycles
56 °C	10 seconds	
72 °C	30 seconds	
↓		
72 °C	2 minutes	


- 3) Check the insert cDNA.
- 4) Culture the transformed *E. coli*.
- 5) Purify the plasmids from transformed *E. coli* and check the sequences using a DNA sequencer.

[Example of use]

The detection of the amplified cDNA by capillary electrophoresis.

Negative Control : mouse IgG

Test sample : Anti Human Ago2 Monoclonal Antibody

Size distribution of synthesized cDNA by using this kit.

The amount and size distribution of cDNA were determined with the Bioanalyzer DNA 1000 Kit. The cDNA were specifically synthesized from anti-Ago2 immunoprecipitates. Cell line; HeLa, Cell numbers ; 1×10^7 cells.

[Reference]

Hayashida, Y., *et al.* : A useful approach to total analysis of RISC-associated RNA., *BMC Research Notes*, 2, 169 (2009)

[Package]

10 reactions.

[Expiration date]

Indicated on the label.

[Storage]

Keep at -20°C .

[Patent]

Target mRNA Cloning Kit *Wako* is patent pending.

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